

How to use

Oralieve 360° Tooth and Mouth Brush

- ✔ 360° bristles to clean multiple surfaces of the mouth simultaneously
- ✔ A small head to get to hard to reach areas of the mouth
- ✔ Soft bristles to provide a gentle clean to sore or sensitive mouths
- ✔ An easy-to-grip handle
- ✔ Can be used to apply Oralieve Moisturising Mouth Gel
- ✔ Holds fluids so can be used for hydration
- ✔ Can be used to clean teeth and remove secretions safely
- ✔ Available in adult and child sizes

Directions for Use

- Move the brush up and down in the mouth like a standard toothbrush.
- The brush can be rotated in the mouth to provide an effective clean on the teeth and gums in hard to reach areas.
- Rinse after use and leave to air dry.

Do not bite or chew the brush. Replace every 3 months or before if the head is showing signs of wear or rotates.

Procurement information

To order through the NHS Supply Chain, use the following information:

Adult 360° Toothbrush: IMC010 PIP Code: 4198115
Child 360° Toothbrush: ILA1100 PIP Code: 4198107

To order directly from Oralieve UK Ltd, email accounts@oralieve.co.uk or phone 01582 439 122.

www.oralieve.co.uk